



QUALITY CONTROL CERTIFICATE

KIT NAME CALPROTECTIN - CHORUS

Code	86600/C03		
Lot No.	949	2024-06	2023-06

REAGENTS

Label	Description	Lot.	Exp.	NO
D9500	CALPROTECTIN	949	2024-06	2X6
D9510	CALPROTECTIN	949	2024-06	1X4
B9510	CALIBRATORE CALPROTECTIN	949	2025-01	1
C9510	CONTROL CALPROTECTIN	949	2025-01	1
G9502	SAMPLE DILUENT	949	2025-06	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	390	≥ 100
Control (mg/kg)	144,1	91,7-208,3

CALIBRATOR

9	0	0	9	4	9	6	2	1	6	2	1	4	7	3	3	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *As...*



CONTROL SERUM

9	0	0	9	4	9	6	1	5	2	8	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *As...*



This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *As...*

DATE 05.09.2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.