



QUALITY CONTROL CERTIFICATE

KIT NAME dsDNA-G - CHORUS

Code	86032		
Lot No.	970	2024-09	2023-09

REAGENTS

Label	Description	Lot.	Exp.	NO
D3281	dsDNA-G	970	2024-09	6X6
B3281	CALIBRATORE dsDNA-G	970	2024-09	1
C3281	CONTROL SERUM dsDNA-G	970	2024-09	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	690	≥ 100
Control Serum (IU/mL)	48,7	31,3-65,4

CALIBRATOR

2	5	1	9	7	0	6	5	1	8	3	1	6	0	3	4	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

251970651831603400
dsDNA-G
CALIBRATOR 0.175 mL [IVD]
LOT 970
 2024-09 2°C - 8°C
DIESSE Monteriggioni Italy CE
B3281

CONTROL SERUM

2	5	1	9	7	0	6	1	2	2	2	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

251970612220000001
dsDNA-G
CONTROL 0.425 mL [IVD]
LOT 970
 2024-09 2°C - 8°C
DIESSE Monteriggioni Italy CE
C3281

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents.....

DATE 20/09/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.