

**QUALITY CONTROL CERTIFICATE**

**KIT NAME** Chorus Chlamydia trachomatis IgG

Code	81261		
Lot No.	970	2024-09	2023-09

**REAGENTS**

Label	Description	Lot.	Exp.	NO
D9640	Chorus Chlamydia trachomatis IgG	970	2024-09	6X6
B9640	CALIBRATOR Chlamydia trachomatis IgG	970	2024-09	1
C9640	POSITIVE CONTROL Chlamydia trachomatis IgG	970	2024-09	1


**CONTROLS**

	Found results	Expected results
Calibrator (mOD)	500	> 100
Positive Control (Index)	1,8	1,2-2,3

**CALIBRATOR**

2	3	3	9	7	0	6	0	0	0	0	0	5	7	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *Abena*



**CONTROL SERUM**

2	3	3	9	7	0	6	1	2	2	3	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *Abena*



*This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *Abena*

DATE *21-09-2023*

**WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!**

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.