

QUALITY CONTROL CERTIFICATE

KIT NAME CYTOMEGALOVIRUS IgM-CHORUS

Code	81011		
Lot No.	962	2024-11	2023-08

REAGENTS

Label	Description	Lot.	Exp.	NO
D2171	Cytomegalovirus IgM	962	2024-11	6X6
B2171	CALIBRATORE Cytomegalovirus IgM	962	2024-11	1
C2171	CONTROL SERUM Cytomegalovirus IgM	962	2024-11	1
G2171	CONFIRMATORY TEST	962	2024-11	5X1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	550	≥ 100
Control Serum (INDEX)	2,1	1,4-2,6
Negative Control (INDEX)	0,0	<0,7

CALIBRATOR

0	7	9	9	6	2	7	0	0	0	0	0	5	1	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

079962700000510000
Cytomegalovirus IgM
CALIBRATOR 0.175 mL [IVD]
LOT 962
 2024-11 2°C-8°C
DIESSSE Monteriggioni Italy CE 0123 B2171

CONTROL SERUM

0	7	9	9	6	2	7	1	4	2	6	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

079962714260000001
Cytomegalovirus IgM
CONTROL+ 0.425 mL [IVD]
LOT 962
 2024-11 2°C-8°C
DIESSSE Monteriggioni Italy CE 0123 C2171

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE 27/09/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.