



**QUALITY CONTROL CERTIFICATE**

**KIT NAME DEAMIDATED GLIADIN-A - CHORUS**

Code	86098		
Lot No.	963	2024-08	2023-08

**REAGENTS**

Label	Description	Lot.	Exp.	NO
D8650	Deamidated Gliadin-A	963	2024-08	6X6
B8650	CALIBRATORE Deamidated Gliadin-A	963	2024-08	1
C8650	CONTROL SERUM Deamidated Gliadin-A	963	2024-08	1

**CONTROLS**

	Found results	Expected results
Calibrator (mOD)	350	≥ 100
Control Serum (AU/mL)	164,1	79,2-400,0

**CALIBRATOR**

2	5	4	9	6	3	6	2	2	4	3	6	2	3	5	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *[Signature]*

254963622436235000  
DEAMIDATED GLIADIN-A  
CALIBRATOR 0.175 mL IVD  
LOT 963  
2024-08 2°C-8°C  
DIESSA Monteriggioni Italy CE

**CONTROL SERUM**

2	5	4	9	6	3	6	3	1	5	7	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *[Signature]*

254963631570000001  
DEAMIDATED GLIADIN-A  
CONTROL 0.425 mL IVD  
LOT 963  
2024-08 2°C-8°C  
DIESSA Monteriggioni Italy CE

*This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *[Signature]*  
DATE 03.10.2023

**WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!**  
This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent,