



QUALITY CONTROL CERTIFICATE

KIT NAME SS-A - CHORUS

Code	86018			
Lot No.	971	2024-12	2023-09	

REAGENTS

Label	Description	Lot.	Exp.	NO
D3181	SS-A	971	2024-12	6X6
B3181	CALIBRATORE SS-A	971	2024-12	1
C3181	CONTROL SERUM SS-A	971	2024-12	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	430	≥ 100
Control Serum (AU/mL)	41,3	23,3-57,0

CALIBRATOR

1	9	3	9	7	1	7	1	1	9	2	9	4	7	4	4	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---


Label checked with barcode reader YES  Operator Signature HB



CONTROL SERUM

1	9	3	9	7	1	7	1	5	2	8	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature HB



*This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... HB

DATE 10.10.2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!  
This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive