



QUALITY CONTROL CERTIFICATE

KIT NAME PR3 - CHORUS

Code	86066		
Lot No.	971	2024-09	2023-09

REAGENTS

Label	Description	Lot.	Exp.	NO
D3512	PR3	971	2024-09	6X6
B3512	CALIBRATORE PR3	971	2024-09	1
C3512	CONTROL SERUM PR3	971	2024-09	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	349	≥ 100
Control Serum (AU/mL)	78,0	56,0-100,0

CALIBRATOR

2	5	2	9	7	1	6	4	1	4	2	9	2	7	4	9	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *HB*

252971641429274900
PR3
CALIBRATOR 0.175 mL IVD
LOT 971
 2024-09 2°C-8°C
DIESE Monteriggioni Italy CE

CONTROL SERUM

2	5	2	9	7	1	6	2	7	4	9	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *HB*

252971627490000001
PR3
CONTROL 0.425 mL IVD
LOT 971
 2024-09 2°C-8°C
DIESE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *HB*

DATE *19.10.2023*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!
This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive