



QUALITY CONTROL CERTIFICATE

KIT NAME CYTOMEGALOVIRUS IgG AVIDITY CHORUS

Code	81092
Lot No.	972
	2024-09
	2023-09

REAGENTS

Label	Description	Lot.	Exp.	NO
D2491	CYTOMEGALOVIRUS IgG AVIDITY	972	2024-09	2x6
C2491	CONTROL SERUM HIGH CMV AVIDITY	972	2024-09	1
G2491	SAMPLE DILUENT	972	2024-09	1

CONTROLS

	Found results	Expected results
Control Serum (%)	90,14%	50-95%



CALIBRATOR


--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--

Label checked with barcode reader YES  Operator Signature .....

CONTROL SERUM

0	8	4	9	7	2	6	5	0	9	5	0	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *AS*



*This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *AS*

DATE 25.10.2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!  
This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.