

QUALITY CONTROL CERTIFICATE

KIT NAME EPSTEIN-BARR EARLY ANTIGEN IgM

Code	81059		
Lot No.	979	2024-10	2023-10

REAGENTS

Label	Description	Lot.	Exp.	NO
D2382	Epstein-Barr Early Antigen IgM	979	2024-10	6x6
B2382	CALIBRATORE Epstein-Barr Early Antigen IgM	979	2024-10	1
C2382	CONTROL SERUM Epstein-Barr Early Antigen IgM	979	2024-10	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	340	≥ 100
Control Serum (Index)	2,9	2,2-4,1
Negative Control (Index)	0,0	<0,7

CALIBRATOR

0	2	9	9	7	9	6	0	0	0	0	0	3	2	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *RB*

02997960000320000
Epstein-Barr Early Antigen IgM
CALIBRATOR 0.178 mL IVD
LOT 979
 2024-10 2°C-8°C
DIESSA Monteriggioni Italy CE

CONTROL SERUM

0	2	9	9	7	9	6	2	2	4	1	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *RB*

02997962241000001
Epstein-Barr Early Antigen IgM
CONTROL 0.425 mL IVD
LOT 979
 2024-10 2°C-8°C
DIESSA Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *Blue*

DATE *09/11/2023*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.