



DIESSSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020

Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME VARICELLA IgG - CHORUS

Code	81078		
Lot No.	974	2024-09	2023-09

REAGENTS

Label	Description	Lot.	Exp.	NO
D2471	VARICELLA IgG	974	2024-09	6x6
B2472	CALIBRATORE VARICELLA IgG	974	2024-09	1
C2471	CONTROL SERUM VARICELLA IgG	974	2024-09	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	428	≥ 100
Control Serum (Index)	2,4	1,9-3,4
Negative Control (Index)	0,0	< 0,7

CALIBRATOR

0	1	4	9	7	4	6	0	0	0	0	0	0	3	8	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

01497460000380000
Varicella IgG
CALIBRATOR 0.425 mL IVD
LOT 974
2024-09 2°C-8°C B2472
DIESSSE Monteriggioni Italy CE

CONTROL SERUM

0	1	4	9	7	4	6	1	9	3	4	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

014974619340000001
Varicella IgG
CONTROL+ 0.425 mL IVD
LOT 974
2024-09 2°C-8°C C2471
DIESSSE Monteriggioni Italy CE

*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE 08/11/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.