



DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020

Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME BORDETELLA PERTUSSIS TOXIN IgA - CHORUS

Code	81274		
Lot No.	964	2024-08	2023-10

REAGENTS

Label	Description	Lot.	Exp	No
D3440	BORDETELLA PERTUSSIS TOXIN IgA	964	2024-08	6X6
B3440	CALIBRATORE BORDETELLA P.TOXIN IgA	964	2024-10	1
C3440	CONTROL SERUM BORDETELLA P.TOXIN IgA	964	2024-10	1


CONTROLS

	Found results	Expected results
CALIBRATOR (mOD)	402	>100
CONTROL SERUM (IU/ml)	44,6	28,1-57,3

CALIBRATOR

9	6	5	9	6	4	6	1	1	8	2	7	4	3	3	8	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *MB*



CONTROL SERUM

9	6	5	9	6	4	6	1	7	3	1	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *MB*



This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *[Signature]*

DATE *07-11-2023*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.