

QUALITY CONTROL CERTIFICATE

KIT NAME HERPES SIMPLEX 1+2 IgM -CHORUS

| | | | |
|---------|-------|---------|---------|
| Code | 81021 | | |
| Lot No. | 979 | 2024-10 | 2023-10 |

REAGENTS

| Label | Description | Lot. | Exp. | NO |
|-------|--------------------------------------|------|---------|-----|
| D2221 | Herpes simplex 1+2 IgM | 979 | 2024-10 | 6x6 |
| B2221 | CALIBRATORE Herpes simplex 1+2 IgM | 979 | 2024-10 | 1 |
| C2221 | CONTROL SERUM Herpes simplex 1+2 IgM | 979 | 2024-10 | 1 |

CONTROLS

| | Found results | Expected results |
|--------------------------|---------------|------------------|
| Calibrator (mOD) | 510 | ≥ 100 |
| Control Serum (INDEX) | 2,5 | 1,8-3,3 |
| Negative Control (INDEX) | 0,1 | <0,7 |

CALIBRATOR

| | | | | | | | | | | | | | | | | | |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| 1 | 5 | 6 | 9 | 7 | 9 | 6 | 0 | 0 | 0 | 0 | 0 | 4 | 0 | 0 | 0 | 0 | 0 |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|


Label checked with barcode reader YES Operator Signature



CONTROL SERUM

| | | | | | | | | | | | | | | | | | |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| 1 | 5 | 6 | 9 | 7 | 9 | 6 | 1 | 8 | 3 | 3 | 0 | 0 | 0 | 0 | 0 | 0 | 1 |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|

Label checked with barcode reader YES Operator Signature



This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE 07.11.2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.