

**QUALITY CONTROL CERTIFICATE**

**KIT NAME** HELICOBACTER PYLORI IgG -CHORUS

<b>Code</b>	81060
<b>Lot No.</b>	981
	2024-10
	2023-10

**REAGENTS**

Label	Description	Lot.	Exp.	NO
D2391	Helicobacter Pylori IgG	981	2024-10	6X6
B2392	CALIBRATORE Helicobacter Pylori IgG	981	2024-10	1
C2391	CONTROL SERUM Helicobacter Pylori IgG	981	2024-10	1

**CONTROLS**

	Found results	Expected results
Calibrator (mOD)	352	≥ 100
Control Serum (UA/mL)	60,0	55,3-141,7
Negative Control (UA/mL)	<5,0	< 5

**CALIBRATOR**

0	0	8	9	8	1	6	1	2	3	4	2	2	8	5	4	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature



**CONTROL SERUM**

0	0	8	9	8	1	6	2	5	4	7	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature



**This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.**

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE 14/11/2023

**WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!**

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.