



QUALITY CONTROL CERTIFICATE

KIT NAME tTg-A - CHORUS

Code	86058		
Lot No.	979	2024-10	2023-10

REAGENTS

Label	Description	Lot.	Exp.	NO
D3431	tTg-A	979	2024-10	6X6
B3431	CALIBRATORE tTg-A	979	2024-10	1
C3431	CONTROL SERUM tTg-A	979	2024-10	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	720	≥ 100
Control Serum (AU/mL)	35,9	21,7-49,1

CALIBRATOR

3	3	9	9	7	9	6	2	1	9	3	0	5	1	3	7	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *HB*

339979621930513700
tTg-A
CALIBRATOR 0,175 mL IVD
LOT 979
2024-10
DIESSE Monteriggioni Italy CE

CONTROL SERUM

3	3	9	9	7	9	6	1	4	2	6	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *HB*

339979614260000001
tTg-A
CONTROL+ 0,425 mL IVD
LOT 979
2024-10
DIESSE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *HB*

DATE *17.11.2023*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.