



QUALITY CONTROL CERTIFICATE

KIT NAME MUMPS IgG

Code	81075		
Lot No.	987	2025-02	2023-11

REAGENTS

Label	Description	Lot.	Exp.	NO
D2451	MUMPS IgG	987	2025-02	6X6
B2452	CALIBRATORE MUMPS IgG	987	2025-02	1
C2451	CONTROL SERUM MUMPS IgG	987	2025-02	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	410	≥ 100
Control Serum (Index)	2,7	2,1-3,9
Negative Control (Index)	0,1	< 0,7

CALIBRATOR

0	1	0	9	8	7	7	0	0	0	0	0	3	3	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

010987700000330000
Mumps IgG
CALIBRATOR 0.425 mL IVD
LOT 987
2025-02 2°C-8°C
DIESSE Monteriggioni Italy CE

CONTROL SERUM

0	1	0	9	8	7	7	2	1	3	9	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

010987721390000001
Mumps IgG
CONTROL+ 0.425 mL IVD
LOT 987
2025-02 2°C-8°C
DIESSE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents.....

DATE 23/11/23

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.