



DIESE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020

Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME MYCOPLASMA PNEUMONIAE IgM - CHORUS

Code	81035
Lot No.	989
	2025-02
	2023-11

REAGENTS

Label	Description	Lot.	Exp.	NO
D2261	MYCOPLASMA PNEUMONIAE IgM	989	2025-02	6x6
B2261	CALIBRATORE MYCOPLASMA PNEUMONIAE IgM	989	2025-02	1
C2261	CONTROL SERUM MYCOPLASMA PNEUMONIAE IgM	989	2025-02	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	608	≥ 100
Control Serum (INDEX)	2,8	1,9-3,6
Negative Control (INDEX)	0,1	≤ 0,7

CALIBRATOR

0	6	7	9	8	9	7	0	0	0	0	0	3	7	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *[Signature]*

06798970000370000
Mycoplasma pneumoniae IgM
CALIBRATOR 0,175 mL [IVD]
LOT 989
2025-02 2°C-8°C
DIESE Monteriggioni Italy CE

CONTROL SERUM

0	6	7	9	8	9	7	1	9	3	6	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *[Signature]*

06798971936000001
Mycoplasma pneumoniae IgM
CONTROL+ 0,425 mL [IVD]
LOT 989
2025-02 2°C-8°C
DIESE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *[Signature]*

DATE 05.12.2013

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.