

QUALITY CONTROL CERTIFICATE

KIT NAME RF-M - CHORUS

Code	86040		
Lot No.	987	2024-11	2023-11

REAGENTS

Label	Description	Lot.	Exp.	NO
D7451	RF-M	987	2024-11	6X6
B7451	CALIBRATORE RF-M	987	2024-11	1
C7451	CONTROL SERUM RF-M	987	2024-11	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	682	≥ 100
Control Serum (IU/mL)	30,5	23,8-115,4

CALIBRATOR

1	6	6	9	8	7	6	1	1	8	2	1	6	0	3	4	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *[Signature]*



CONTROL SERUM

1	6	6	9	8	7	6	1	2	2	2	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *[Signature]*



This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *[Signature]*

DATE 07-12-2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.