



QUALITY CONTROL CERTIFICATE

KIT NAME DEAMIDATED GLIADIN-A - CHORUS

Code	86098		
Lot No.	988	2024-11	2023-11

REAGENTS

Label	Description	Lot.	Exp.	NO
D8650	Deamidated Gliadin-A	988	2024-11	6X6
B8650	CALIBRATORE Deamidated Gliadin-A	988	2024-11	1
C8650	CONTROL SERUM Deamidated Gliadin-A	988	2024-11	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	315	≥ 100
Control Serum (AU/mL)	91,3	45,9-200,0

CALIBRATOR

2	5	4	9	8	8	6	3	2	7	3	9	2	8	5	5	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *Abu*

254988632739285500
 DEAMIDATED GLIADIN-A
CALIBRATOR 0.175 mL IVD
 LOT 988
 2024-11 2°C - 8°C
 DIESSSE Monteriggioni Italy CE

CONTROL SERUM

2	5	4	9	8	8	6	2	5	4	7	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *Abu*

254988625470000001
 DEAMIDATED GLIADIN-A
CONTROL+ 0.425 mL IVD
 LOT 988
 2024-11 2°C - 8°C
 DIESSSE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *Abu*

DATE 18 Oct 2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.