



QUALITY CONTROL CERTIFICATE

KIT NAME MUMPS IgM CHORUS

Code 81076/12

Lot No. 001

2025-01

2024-01

REAGENTS

Label	Description	Lot.	Exp.	NO
D2461	MUMPS IgM	001	2025-01	2x6
B2461	CALIBRATORE MUMPS IgM	001	2025-01	1
C2461	CONTROL SERUM MUMPS IgM	001	2025-01	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	430	≥ 100
Control Serum (Index)	2,0	1,4-2,7
Negative Control (Index)	0,1	< 0,7

13200160000490000
Mumps IgM
CALIBRATOR 0.175 mL IVD
LOT 001
 2025-01 2024-01
DIESSA Monteriggioni Italy CE

CALIBRATOR

1 3 2 0 0 1 6 0 0 0 0 0 4 9 0 0 0 0

Label checked with barcode reader

YES

Operator Signature

CONTROL SERUM

1 3 2 0 0 1 6 1 4 2 7 0 0 0 0 0 0 1

Label checked with barcode reader

YES

Operator Signature

13200161427000001
Mumps IgM
CONTROL+ 0.425 mL IVD
LOT 001
 2025-01 2024-01
DIESSA Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE 07.02.2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.