



QUALITY CONTROL CERTIFICATE

KIT NAME LH

Code	86512		
Lot No.	990	2024-11	2024-01

REAGENTS

Label	Description	Lot.	Exp.	NO
D3250	LH	990	2024-11	6X6
B3250	CALIBRATOR LH	990	2024-11	1
C3250	CONTROL SERUM LH	990	2024-11	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	4,30	≥ 100
Control Serum (mIU/ml)	21,6	12,8-29,2

CALIBRATOR

9	2	2	9	9	0	6	7	1	4	3	0	2	7	4	9	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *Olivero*

922990671430274900
LH
CALIBRATOR 0.500 mL IVD
LOT 990
2024-11
2°C-8°C
DIESE Monteriggioni Italy CE

CONTROL SERUM

9	2	2	9	9	0	6	2	6	4	8	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *Olivero*

922990626480000001
LH
CONTROL+ 1 mL IVD
LOT 990
2024-11
2°C-8°C
DIESE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *Olivero*

DATE 08-02-2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!
This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.