

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME POLIOVIRUS IgG - CHORUS

Code	81264		
Lot No.	987	2024-11	2024-01

REAGENTS

Label	Description	Lot.	Exp.	NO
D9630	Poliovirus IgG	987	2024-11	6X6
B9630	CALIBRATORE Poliovirus IgG	987	2024-12	1
C9630	CONTROL SERUM Poliovirus IgG	987	2024-12	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	473	≥ 100
Control (Index)	1,4	1,1-2,1

CALIBRATOR

9	0	1	9	8	7	6	0	0	0	0	0	6	2	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader **YES** Operator Signature SS

90198760000620000
Poliovirus IgG
CALIBRATOR 0.176 mL [IVD]
LOT 987
2024-12 2°C-8°C
DIESSE Monteriggioni Italy CE

CONTROL SERUM

9	0	1	9	8	7	6	1	1	2	1	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader **YES** Operator Signature SS

90198761121000001
Poliovirus IgG
CONTROL+ 0.425 mL [IVD]
LOT 987
2024-12 2°C-8°C
DIESSE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents SS

DATE 16.02.24

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.