

QUALITY CONTROL CERTIFICATE

KIT NAME Epstein-Barr VCA IgM II - CHORUS

Code	81054		
Lot No.	001	2025-04	2024-02

REAGENTS

Label	Description	Lot.	Exp.	NO
D9490	Epstein-Barr VCA IgM II	001	2025-04	6x6
B9490	CALIBRATOR Epstein-Barr VCA IgM II	001	2025-05	1
C9490	CONTROL SERUM Epstein-Barr VCA IgM II	001	2025-05	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	450	≥ 100
Control Serum (INDEX)	3,4	2,2-4,1
Negative Control (INDEX)	0,1	< 0,7

CALIBRATOR

0	3	0	0	0	0	1	7	0	0	0	0	0	0	3	2	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature SB

03000170000320000
Epstein-Barr VCA IgM II
CALIBRATOR 0.175 mL IVD
LOT 001
2025-05 2°C-8°C
DIESSE Monteriggioni Italy CE

CONTROL SERUM

0	3	0	0	0	0	1	7	2	2	4	1	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature SB

03000172241000001
Epstein-Barr VCA IgM II
CONTROL+ 0.425 mL IVD
LOT 001
2025-05 2°C-8°C
DIESSE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents: SB

DATE 10.02.2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.