

**QUALITY CONTROL CERTIFICATE**

**KIT NAME** Epstein-Barr EBNA IgG -CHORUS

Code	81057		
Lot No.	002	2025-01	2024-01

**REAGENTS**

Label	Description	Lot.	Exp.	NO
D2361	Epstein-Barr EBNA IgG	002	2025-01	6X6
B2362	CALIBRATORE Epstein-Barr EBNA IgG	002	2025-01	1
C2361	CONTROL SERUM Epstein-Barr EBNA IgG	002	2025-01	1

**CONTROLS**

	Found results	Expected results
Calibrator (mOD)	400	≥ 100
Control Serum (INDEX)	2,9	2,0-3,7
Negative Control (INDEX)	0,0	<0,7

**CALIBRATOR**

0	0	6	0	0	2	6	0	0	0	0	0	3	5	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature MB

00600260000350000  
Epstein-Barr EBNA IgG  
CALIBRATOR 0.425 mL IVD  
LOT 002  
2025-01 2°C 8°C  
DIESE Monteriggioni Italy CE

**CONTROL SERUM**

0	0	6	0	0	2	6	2	0	3	7	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature MB

006002620370000001  
Epstein-Barr EBNA IgG  
CONTROL+ 0.425 mL IVD  
LOT 002  
2025-01 2°C 8°C  
DIESE Monteriggioni Italy CE

**This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.**

APPROVED RESPONSIBLE HEAD, Quality Control Reagents [Signature]

DATE 20/02/2024

**WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!**

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.