



DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020

Sede Legale  
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME TOXOCARA IgG -CHORUS

Code	81280		
Lot No.	016	2025-03	2024-03

REAGENTS

Label	Description	Lot.	Exp.	NO
D0020	TOXOCARA IgG	016	2025-03	2X6
B0020	CALIBRATORE Toxocara IgG	016	2025-03	1
C0020	CONTROL SERUM Toxocara IgG	016	2025-03	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	540	≥ 100
Control (Index)	1,9	1,2-2,3

CALIBRATOR

9	3	3	0	1	6	6	0	0	0	0	0	5	8	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature HB

933016600000580000  
Toxocara IgG  
CALIBRATOR 0,300 mL IVD  
LOT 016  
2025-03 2°C-8°C  
DIESSE Monteriggioni Italy CE

CONTROL SERUM

9	3	3	0	1	6	6	1	2	2	3	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature HB

933016612230000001  
Toxocara IgG  
CONTROL 0,600 mL IVD  
LOT 016  
2025-03 2°C-8°C  
DIESSE Monteriggioni Italy CE

**This batch has been manufactured and checked according to quality requirements.  
Hereby, it is released to distribution.**

APPROVED RESPONSIBLE HEAD, Quality Control Reagents [Signature]

DATE 13.03.2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.