



QUALITY CONTROL CERTIFICATE

KIT NAME MUMPS IgG

Code	81075			
Lot No.	988	2025-02	2023-11	

REAGENTS

Label	Description	Lot.	Exp.	NO
D2451	MUMPS IgG	988	2025-02	6X6
B2452	CALIBRATORE MUMPS IgG	988	2025-02	1
C2451	CONTROL SERUM MUMPS IgG	988	2025-02	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	360	≥ 100
Control Serum (Index)	3,5	2,1-4,0
Negative Control (Index)	0,0	< 0,7

CALIBRATOR

0	1	0	9	8	8	7	0	0	0	0	0	3	3	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature MB

010988700000330000  
Mumps IgG  
CALIBRATOR 0.425 mL [IVD]  
LOT 988  
2025-02 2°C 8°C  
DIESSSE Monteriggioni Italy CE

CONTROL SERUM

0	1	0	9	8	8	7	2	1	4	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature MB

010988721400000001  
Mumps IgG  
CONTROL+ 0.425 mL [IVD]  
LOT 988  
2025-02 2°C 8°C  
DIESSSE Monteriggioni Italy CE

**This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.**

APPROVED RESPONSIBLE HEAD, Quality Control Reagents [Signature]

DATE 13/03/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.