

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME MEASLES IgM CHORUS

Code	81073		
Lot No.	979	2025-01	2023-10

REAGENTS

Label	Description	Lot.	Exp.	NO
D2442	MEASLES IgM	979	2025-01	6x6
B2442	CALIBRATORE MEASLES IgM	979	2025-01	2
C2442	CONTROL SERUM MEASLES IgM	979	2025-01	2

CONTROLS

	Found results	Expected results
Calibrator (mOD)	570	≥ 100
Control Serum (Index)	1,5	1,3-2,4
Negative Control (Index)	0,0	< 0,7

CALIBRATOR

1	3	3	9	7	9	7	0	0	0	0	0	5	5	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *[Signature]*

13397970000550000
Measles IgM
CALIBRATOR 0.175 mL IVD
LOT 979
 2025-01 2°C-8°C
DIESSE Monteriggioni Italy CE

CONTROL SERUM

1	3	3	9	7	9	7	1	3	2	4	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *[Signature]*

133979713240000001
Measles IgM
CONTROL+ 0.425 mL IVD
LOT 979
 2025-01 2°C-8°C
DIESSE Monteriggioni Italy CE

*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *[Signature]*

DATE 11.03.2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.