



DIESSSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020

Sede Legale  
20144 Milano (Italy) Via A. Solari 19 sc. 6

**QUALITY CONTROL CERTIFICATE**

**KIT NAME** CHORUS PHOSPHOLIPID SCREEN IgG

Code	86106		
Lot No.	016	2025-03	2024-03

**REAGENTS**

Label	Description	Lot.	Exp.	NO
D4010	Phospholipid Screen IgG	016	2025-03	6X6
B4010	CALIBRATORE Phospholipid Screen IgG	016	2025-03	1
C4010	CONTROL SERUM Phospholipid Screen IgG	016	2025-03	1

**CONTROLS**

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	540	≥ 100
Control (AU/ml)	33,8	25,0-55,2

**CALIBRATOR**

9	7	4	0	1	6	6	6	1	4	2	2	4	4	4	3	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature MB

974016661422444300  
Phospholipid Screen-G  
CALIBRATOR 0.450 mL IVD  
LOT 016  
 2025-03 2°C-8°C  
DIESSSE Monteriggioni Italy CE  
B4010

**CONTROL SERUM**

9	7	4	0	1	6	6	1	6	3	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature MB

974016616300000001  
Phospholipid Screen-G  
CONTROL+ 0.800 mL IVD  
LOT 016  
 2025-03 2°C-8°C  
DIESSSE Monteriggioni Italy CE  
C4010

*This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE 20/03/2024

**WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!**

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.