



**QUALITY CONTROL CERTIFICATE**

**KIT NAME** ANA-Screen - CHORUS

Code	86014			
Lot No.	016	2025-03	2024-03	

**REAGENTS**

Label	Description	Lot.	Exp.	NO
D3141	ANA-Screen	016	2025-03	6X6
B3141	CALIBRATORE ANA-Screen	016	2025-03	1
C3141	CONTROL SERUM ANA-Screen	016	2025-03	1

**CONTROLS**

	Found results	Expected results
Calibrator (mOD)	360	≥ 100
Control Serum (Ndx)	3,4	2,3-4,3

**CALIBRATOR**

1	1	5	0	1	6	6	0	0	0	0	0	3	0	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature TS

11501660000300000  
ANA-Screen  
CALIBRATOR 0.176 mL IVD  
LOT 016  
2025-03  
DIESSE Monteriggioni Italy CE

**CONTROL SERUM**

1	1	5	0	1	6	6	2	3	4	3	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature TS

115016623430000001  
ANA-SCREEN  
CONTROL+ 0.425 mL IVD  
LOT 016  
2025-03  
DIESSE Monteriggioni Italy CE

*This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents. [Signature]

DATE 25/03/2024

**WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!**

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.