

QUALITY CONTROL CERTIFICATE

KIT NAME INFLUENZA A IgG -CHORUS

Code	81190		
Lot No.	016	2025-03	2024-03

REAGENTS

Label	Description	Lot.	Exp.	NO
D7511	INFLUENZA A IgG	016	2025-03	6X6
B7511	CALIBRATORE INFLUENZA A IgG	016	2025-03	1
C7511	CONTROL SERUM INFLUENZA A IgG	016	2025-03	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	410	≥ 100
Control Serum (Index)	3,0	1,9-3,5
Negative Control (Index)	0,0	≤ 0,7

CALIBRATOR

0	7	0	0	1	6	6	0	0	0	0	0	3	7	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature MB

07001660000370000
Influenza A IgG
CALIBRATOR 0.175 mL IVD
LOT 016
 2025-03 2°C-8°C
DIESSE Monteriggioni Italy CE

CONTROL SERUM

0	7	0	0	1	6	6	1	9	3	5	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature MB

0700166193500000001
Influenza A IgG
CONTROL+ 0.425 mL IVD
LOT 016
 2025-03 2°C-8°C
DIESSE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents MB

DATE 29.03.2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.