



QUALITY CONTROL CERTIFICATE

KIT NAME CARDIOLIPIN-G - CHORUS

Code	86046
Lot No.	025-A 2025-04 2024-04

REAGENTS

Label	Description	Lot.	Exp.	NO
D3311	CARDIOLIPIN-G	025	2025-04	6X6
B3311	CALIBRATORE CARDIOLIPIN-G	025-1	2025-04	1
C3311	CONTROL SERUM CARDIOLIPIN-G	025-1	2025-04	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	430	≥ 100
Control Serum (GPLU/mL)	30,8	22,5-43,1

CALIBRATOR

3	3	5	0	2	5	6	3	2	2	3	8	4	5	5	9	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

335025632238455900
Cardiolipin-G
CALIBRATOR 0.175 mL IVD
LOT 025-1
2025-04 2°C - 8°C
DIESSE Monteriggioni Italy CE

CONTROL SERUM

3	3	5	0	2	5	6	1	6	2	9	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

335025616290000001
Cardiolipin-G
CONTROL+ 0.425 mL IVD
LOT 025-1
2025-04 2°C - 8°C
DIESSE Monteriggioni Italy CE

*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE 08-05-2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!
This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.