

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME CALPROTECTIN 3RD GENERATION- CHORUS

Code	86603		
Lot No.	017	2025-03	2024-03

REAGENTS

Label	Description	Lot.	Exp.	NO
DB430	CALPROTECTIN 3 RD GENERATION	017	2025-03	6X6
B8430	CALIBRATORE CALPROTECTIN 3 RD GENERATION	017	2025-03	1
C8430	CONTROL CALPROTECTIN 3 RD GENERATION	017	2025-03	1
G8430	SAMPLE DILUENT/NEGATIVE CONTROL CALPROTECTIN 3 RD GENERATION	017	2025-09	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	350	≥ 100
Control (mg/Kg)	564,6	309,2-800,0
sample dil/neg control Calp 3 RD (mg/Kg)	24,0	< 51,0

CALIBRATOR

9	9	3	0	1	7	6	3	2	3	3	6	3	3	4	5	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---


Label checked with barcode reader YES ☒ Operator Signature



CONTROL SERUM

9	9	3	0	1	7	6	2	1	3	9	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES ☒ Operator Signature



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED BY RESPONSIBLE HEAD, Quality Control Reagents.

DATE 09.05.2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.