



DIESSSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020

Sede Legale  
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME MEASLES IgG CHORUS

Code	81072
Lot No.	025
	2025-07
	2024-04

REAGENTS

Label	Description	Lot.	Exp.	NO
D2431	MEASLES IgG	025	2025-07	6X6
B2432	CALIBRATORE MEASLES IgG	025	2025-07	1
C2431	CONTROL SERUM MEASLES IgG	025	2025-07	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	380	≥ 100
Control Serum (Index)	3,2	2,3-4,2
Negative Control (Index)	0,0	< 0,7

CALIBRATOR

0	0	9	0	2	5	7	0	0	0	0	0	3	1	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *Così*

009025700000310000  
Measles IgG  
CALIBRATOR 0.425 mL IVD  
LOT 025  
2025-07  
DIESSSE Monteriggioni Italy CE

CONTROL SERUM

0	0	9	0	2	5	7	2	3	4	2	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *Così*

009025723420000001  
Measles IgG  
CONTROL 0.425 mL IVD  
LOT 025  
2025-07  
DIESSSE Monteriggioni Italy CE

**This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.**

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *Così*

DATE 13.05.2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.