

DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554



QUALITY CONTROL CERTIFICATE

KIT NAME HELICOBACTER PYLORI IgG -CHORUS

Code 81060  
Lot No. 025 2025-04 2024-04

REAGENTS

Label	Description	Lot.	Exp.	NO
D2392	Helicobacter Pylori IgG	025	2025-04	6X6
B2393	CALIBRATORE Helicobacter Pylori IgG	025	2025-04	1
C2392	CONTROL SERUM Helicobacter Pylori IgG	025	2025-04	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	348	≥ 100
Control Serum (UA/mL)	103,5	53,1-110,7
Negative Control (UA/mL)	<5,0	< 5

CALIBRATOR

0 0 8 0 2 5 6 2 2 0 3 6 3 4 6 4 0 0

Label checked with barcode reader YES  Operator Signature *[Signature]*

008025622036346400  
HELICOBACTER PYLORI IgG  
CALIBRATOR 0.425 mL IVD  
LOT 025  
2025-04 2°C 4°C  
DIESSE Monteriggioni Italy CE 0123 B2393

CONTROL SERUM

0 0 8 0 2 5 6 2 1 3 9 0 0 0 0 0 0 1

Label checked with barcode reader YES  Operator Signature *[Signature]*

008025621390000001  
HELICOBACTER PYLORI IgG  
CONTROL+ 0.425 mL IVD  
LOT 025  
2025-04 2°C 4°C  
DIESSE Monteriggioni Italy CE 0123 C2392

**This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.**

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *[Signature]*

DATE 16-05-2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.