

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME **HELICOBACTER PYLORI Ag - CHORUS**

Code **81063**

Lot No. **027**

2025-04

2024-05

REAGENTS

Label	Description	Lot	Exp.	NO
D6960	Helicobacter pylori Ag	027	2025-04	2X6
B6960	CALIBRATORE Helicobacter pylori Ag	027	2025-08	2
C6960	CONTROL SERUM Helicobacter pylori Ag	027	2025-08	3
G6960	H. pylori Ag Stool Diluent/Negative Control	E0330200	2025-09	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	299	≥ 100
Control + (Index)	3,9	2,7-5,1
Stool Diluent/Negative Control (Index)	0,2	$< 0,5$


CALIBRATOR

9 3 2 0 2 7 6 0 0 0 0 0 2 6 0 0 0 0

Label checked with barcode reader

YES ☒

Operator Signature




CONTROL SERUM

9 3 2 0 2 7 6 2 7 5 1 0 0 0 0 0 0 1

Label checked with barcode reader

YES ☒

Operator Signature



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE

26.05.2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.