

DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554



QUALITY CONTROL CERTIFICATE

KIT NAME INFLUENZA A IgA -CHORUS

Code	81192		
Lot No.	033	2025-05	2024-05

REAGENTS

Label	Description	Lot.	Exp.	NO
D7522	INFLUENZA A IgA	033	2025-05	6X6
B7522	CALIBRATORE INFLUENZA A IgA	033	2025-05	1
C7522	CONTROL SERUM INFLUENZA A IgA	033	2025-05	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	500	≥ 100
Control Serum (Ndx)	2,2	1,4-2,6
Negative Control (Ndx)	0,3	≤ 0,7

CALIBRATOR

0	6	9	0	3	3	6	0	0	0	0	0	5	0	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature MB

06903360000500000  
INFLUENZA A IgA  
CALIBRATOR 0.178 mL IVD  
LOT 033  
2025-05 2°C-8°C  
DIESSE Monteriggioni Italy CE 0123

CONTROL SERUM

0	6	9	0	3	3	6	1	4	2	6	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature MB

069033614260000001  
INFLUENZA A IgA  
CONTROL SERUM 0.425 mL IVD  
LOT 033  
2025-05 2°C-8°C  
DIESSE Monteriggioni Italy CE 0123

*This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents [Signature]  
DATE 28/05/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!  
This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.