



QUALITY CONTROL CERTIFICATE

KIT NAME **Epstein-Barr VCA IgM II - CHORUS**

Code	81054		
Lot No.	005	2025-04	2024-02

REAGENTS

Label	Description	Lot	Exp.	NO
D9490	Epstein-Barr VCA IgM II	005	2025-04	6x6
B9490	CALIBRATOR Epstein-Barr VCA IgM II	005	2025-05	1
C9490	CONTROL SERUM Epstein-Barr VCA IgM II	005	2025-05	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	490	≥ 100
Control Serum (INDEX)	2,6	1,7-3,1
Negative Control (INDEX)	0,1	< 0,7

CALIBRATOR

0 3 0 0 0 5 7 0 0 0 0 0 4 2 0 0 0 0

Label checked with barcode reader YES Operator Signature

03000570000420000
Epstein-Barr VCA IgM II
CALIBRATOR 0.175 mL IVD
LOT 005
2025-05 2°C-8°C BB490
DIESSE Monteriggioni Italy CE

CONTROL SERUM

0 3 0 0 0 5 7 1 7 3 1 0 0 0 0 0 0 0 1

Label checked with barcode reader YES Operator Signature

030005717310000001
Epstein-Barr VCA IgM II
CONTROL 0.425 mL IVD
LOT 005
2025-05 2°C-8°C C9490
DIESSE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents.....

DATE **30-05-2024**

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.