



DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020

Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME HSV 1 SCREEN - CHORUS

Code 81018

Lot No. 033

2025-08

2024-05

REAGENTS

Label	Description	Lot	Exp.	NO
D2191	HSV 1 SCREEN	033	2025-08	6x6
B2191	CALIBRATORE HSV 1 SCREEN	033	2025-09	1
C2191	CONTROL SERUM HSV 1 SCREEN	033	2025-08	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

CALIBRATOR (mOD)	520	≥ 100
CONTROL SERUM (Index)	3,0	1,3-8,0
NEGATIVE CONTROL (Index)	0,5	< 0,7

CALIBRATOR

0 3 3 0 3 3 7 0 0 0 0 0 5 8 0 0 0 0


Label checked with barcode reader YES Operator Signature *Cher...*



CONTROL SERUM

0 3 3 0 3 3 7 1 3 8 0 0 0 0 0 0 0 1

Label checked with barcode reader YES Operator Signature *Cher...*



This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *Cher...*

DATE 06.06.2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.