

QUALITY CONTROL CERTIFICATE

KIT NAME VARICELLA IgM CHORUS

Code	81079		
Lot No.	033	2025-08	2024-05

REAGENTS

Label	Description	Lot.	Exp.	NO
D2481	VARICELLA IgM	033	2025-08	6x6
B2481	CALIBRATORE VARICELLA IgM	033	2025-08	1
C2481	CONTROL SERUM VARICELLA IgM	033	2025-08	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	430	≥ 100
Control Serum (Index)	3,3	2,1-4,0
Negative Control (Index)	0,0	< 0,7

CALIBRATOR

1	3	4	0	3	3	7	0	0	0	0	0	3	3	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature

134033700000330000  
Varicella IgM  
CALIBRATOR 0.175 mL IVD  
LOT 033  
2025-08  
DIESSA Monteriggioni Italy CE

CONTROL SERUM

1	3	4	0	3	3	7	2	1	4	0	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature

134033721400000001  
Varicella IgM  
CONTROL 0.425 mL IVD  
LOT 033  
2025-08  
DIESSA Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements.  
Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE 03/06/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.