



DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020

Sede Legale  
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME EPSTEIN-BARR EARLY ANTIGEN IgG

Code	81058		
Lot No.	033	2025-08	2024-05

REAGENTS

Label	Description	Lot	Exp.	NO
D2372	Epstein-Barr Early Antigen IgG	033	2025-08	6x6
B2373	CALIBRATORE Epstein-Barr Early Antigen IgG	033	2025-08	1
C2372	CONTROL SERUM Epstein-Barr Early Antigen IgG	033	2025-08	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	345	≥ 100
Control Serum (Index)	3,0	2,3-4,4
Negative Control (Index)	0,0	< 0,7

CALIBRATOR

0	0	7	0	3	3	7	0	0	0	0	0	0	3	0	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature ..... MB

007033700000300000  
Epstein-Barr Early Antigen IgG  
CALIBRATOR 0.425 mL IVD  
LOT 033  
2025-08 2°C 8°C  
DIESSE Monteriggioni Italy CE

CONTROL SERUM

0	0	7	0	3	3	7	2	3	4	4	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature ..... MB

007033723440000001  
Epstein-Barr Early Antigen IgG  
CONTROL 0.425 mL IVD  
LOT 033  
2025-08 2°C 8°C  
DIESSE Monteriggioni Italy CE

*This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents.....

DATE 05/06/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.