

DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554



## QUALITY CONTROL CERTIFICATE

KIT NAME DEAMIDATED GLIADIN-G - CHORUS

Code 86100

Lot No. 026-A

2025-04

2024-04

### REAGENTS

| Label | Description                        | Lot   | Exp.    | NO  |
|-------|------------------------------------|-------|---------|-----|
| D8661 | Deamidated Gliadin-G               | 026   | 2025-04 | 6X6 |
| B8661 | CALIBRATORE Deamidated Gliadin-G   | 026-1 | 2025-04 | 1   |
| C8661 | CONTROL SERUM Deamidated Gliadin-G | 026-1 | 2025-04 | 1   |


### CONTROLS

|                       | Found results | Expected results |
|-----------------------|---------------|------------------|
| Calibrator (mOD)      | 530           | $\geq 100$       |
| Control Serum (AU/mL) | 44,7          | 27,7-72,7        |

### CALIBRATOR

2 5 5 0 2 6 6 1 2 1 3 2 5 0 5 4 0 0

Label checked with barcode reader YES ☒ Operator Signature *Green*



### CONTROL SERUM

2 5 5 0 2 6 6 1 4 2 6 0 0 0 0 0 0 0 1

Label checked with barcode reader YES ☒ Operator Signature *Green*



*This batch has been manufactured and checked according to quality requirements.  
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *Green*

DATE *10.06.2024*

### WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.