

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554



QUALITY CONTROL CERTIFICATE

KIT NAME Epstein-Barr VCA IgM II - CHORUS

Code 81054

Lot No. 007

2025-05

2024-05

REAGENTS

Label	Description	Lot.	Exp.	NO
D9490	Epstein-Barr VCA IgM II	007	2025-05	6x6
B9490	CALIBRATOR Epstein-Barr VCA IgM II	007	2025-05	1
C9490	CONTROL SERUM Epstein-Barr VCA IgM II	007	2025-05	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	470	≥ 100
Control Serum (INDEX)	2,3	1,6-2,9
Negative Control (INDEX)	0,1	$< 0,7$

CALIBRATOR

0 3 0 0 0 7 7 0 0 0 0 0 4 5 0 0 0 0

Label checked with barcode reader

YES ☒

Operator Signature *MB*



CONTROL SERUM

0 3 0 0 0 7 7 1 6 2 9 0 0 0 0 0 0 0 1

Label checked with barcode reader

YES ☒

Operator Signature *MB*



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents. *Abun*

DATE 12/06/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.