

DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554



QUALITY CONTROL CERTIFICATE

KIT NAME CARDIOLIPIN-G - CHORUS

Code	86046		
Lot No.	034	2025-05	2024-05

REAGENTS

Label	Description	Lot	Exp.	NO
D3312	CARDIOLIPIN-G	034	2025-05	6X6
B3312	CALIBRATORE CARDIOLIPIN-G	034	2025-05	1
C3312	CONTROL SERUM CARDIOLIPIN-G	034	2025-05	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	585	≥ 100
Control Serum (GPLU/mL)	22,5	19,0-32,0

CALIBRATOR

3	3	5	0	3	4	6	4	2	5	4	0	5	1	6	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature

335034642540516000  
CARDIOLIPIN-G  
CALIBRATOR 0.175 mL IVD  
LOT 034  
2025-05  
DIESSE Monteriggioni Italy  
CE 0123  
E3312

CONTROL SERUM

3	3	5	0	3	4	6	1	4	2	6	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature

335034614260000001  
CARDIOLIPIN-G  
CONTROL 0.425 mL IVD  
LOT 034  
2025-05  
DIESSE Monteriggioni Italy  
CE 0123  
C3312

*This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents.

DATE 14/06/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.