



DIESE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554

QUALITY CONTROL CERTIFICATE

KIT NAME HERPES SIMPLEX 2 IgG RECOMBINANT - CHORUS

Code	81029		
Lot No.	026	2025-04	2024-04

REAGENTS

Label	Description	Lot	Exp.	NO
D8690	Herpes simplex 2 IgG Recombinant	026	2025-04	6X6
B8690	CALIBRATORE Herpes simplex 2 IgG Recombinant	026	2025-04	1
C8690	CONTROL SERUM Herpes simplex 2 IgG Recombinant	026	2025-04	1


CONTROLS

	Found results	Expected results
Calibrator (mOD)	540	≥ 100
Control Serum (INDEX)	1,7	1,2-2,1
Negative Control (INDEX)	0,0	<0,7

CALIBRATOR

2	6	0	0	2	6	6	0	0	0	0	0	6	1	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature MIB



CONTROL SERUM

2	6	0	0	2	6	6	1	2	2	1	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature MIB



This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents [Signature]

DATE 20/06/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.