

DIESSSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554



QUALITY CONTROL CERTIFICATE

KIT NAME MYCOPLASMA PNEUMONIAE IgM - CHORUS

Code	81035		
Lot No.	013	2025-05	2024-02

REAGENTS

Label	Description	Lot.	Exp.	NO
D2262	MYCOPLASMA PNEUMONIAE IgM	013	2025-05	6x6
B2262	CALIBRATORE MYCOPLASMA PNEUMONIAE IgM	013	2025-06	1
C2262	CONTROL SERUM MYCOPLASMA PNEUMONIAE IgM	013	2025-06	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	590	≥ 100
Control Serum (INDEX)	1,4	0,9-1,7
Negative Control (INDEX)	0,1	≤ 0,7

CALIBRATOR

0	6	7	0	1	3	7	0	0	0	0	0	8	0	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *MB*



CONTROL SERUM

0	6	7	0	1	3	7	0	9	1	7	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *MB*



This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *[Signature]*

DATE *21.06.2024*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!
This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.