

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554



QUALITY CONTROL CERTIFICATE

KIT NAME RUBELLA IgG - CHORUS

Code 81030

Lot No. 034

2025-05

2024-05

REAGENTS

Label	Description	Lot.	Exp.	NO
D2231	RUBELLA IgG	034	2025-05	6X6
B2232	CALIBRATORE RUBELLA IgG	034	2025-05	1
C2231	CONTROL SERUM RUBELLA IgG	034	2025-05	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	460	≥ 100
Control Serum (IU/mL)	27,2	18,3-41,7
Negative Control (IU/mL)	$<5,0$	< 5

CALIBRATOR

0	8	1	0	3	4	6	5	1	9	3	1	4	8	5	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES ☒

Operator Signature MB



CONTROL SERUM

0	8	1	0	3	4	6	1	5	2	7	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES ☒

Operator Signature MB



This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE

24/06/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.