

QUALITY CONTROL CERTIFICATE

KIT NAME ANTI-CCP - CHORUS

Code	86094		
Lot No.	042	 2025-06	 2024-06

REAGENTS

Label	Description	Lot	Exp.	NO
D8671	ANTI-CCP	042	2025-06	6x6
B8671	CALIBRATORE ANTI-CCP	042	2025-06	1
C8671	CONTROL SERUM ANTI-CCP	042	2025-06	1

CONTROLS


	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	522	≥ 100
Control Serum (AU/mL)	41,8	30,0-63,6

CALIBRATOR

2	5	6	0	4	2	6	2	2	4	4	6	4	3	6	6	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES ☒ Operator Signature HB



CONTROL SERUM

2	5	6	0	4	2	6	1	6	3	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES ☒ Operator Signature HB



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents [Signature]

DATE 25/06/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.