

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554



QUALITY CONTROL CERTIFICATE

KIT NAME PR3 - CHORUS

Code	86066		
Lot No.	042	2025-06	2024-06

REAGENTS

Label	Description	Lot	Exp.	NO
D3512	PR3	042	2025-06	6X6
B3512	CALIBRATORE PR3	042	2025-06	1
C3512	CONTROL SERUM PR3	042	2025-06	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	410	≥ 100
Control Serum (AU/mL)	57,5	39,4-69,5

CALIBRATOR

2	5	2	0	4	2	6	6	1	4	3	0	3	7	5	1	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature MB

252042661430375100
PR3
CALIBRATOR 0.175 mL IVD
LOT 042
2025-06 2°C-8°C
DIESSE Monteriggioni Italy CE

CONTROL SERUM

2	5	2	0	4	2	6	1	9	3	5	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature MB

252042619350000001
PR3
CONTROL 0.425 mL IVD
LOT 042
2025-06 2°C-8°C
DIESSE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents Aburo

DATE 28/06/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.